SIMPLE

SIMPLE Hip Brace N3

REF 21003

Accessories list
Instructions ----- 1 copy
Adjustment Tool ----- 1 piece

SIMPLE Hüftorthese N3
SIMPLE – orthèse de hanche N3
SIMPLE heuporthese N3
Ortesi d'anca SIMPLE N3
Ortesis de cadera SIMPLE N3
Órtese de quadril SIMPLE N3
Orteza biodra SIMPLE N3
Kyčelní ortéza SIMPLE N3
Bedrová ortéza SIMPLE N3
Ортез тазобедренного сустава SIMPLE N3
SIMPLE 고관절 보조기 N3
简单髋部固定器

The detailed symbol glossary and The product is not allowed to be reused can be found on
<https://www.simplemedical.cn/support-detail/148.htm>

[Responsible Person in the EU]

Company: MedNet EC-REP GmbH
Address: Borkstrasse 10,
48163 Münster, Germany
E-Mail: contact@mednet-ecrep.com

Address: 909 ChuangXiang Building 2,
Minzhi Street, Longhua District,
Shenzhen, China. 518131
Tel.: +86-755-86035863
Web: www.simplemedical.cn
E-mail: service@simplegroup.net

21003 V0.3/2026.05

*Please refer to the instructions before usage.

English

Thank you for purchasing this product.
Please read these instructions carefully before use.
Keep this manual for future reference.

Product Overview

Product Description: This product is made of aluminum alloy, stainless steel, sponge, and plastic. It is worn on the outside of the body. It supports and stabilizes the hip by limiting movement.
Note: This product does not correct bone or joint conditions.

Reprocessing instructions: The product features an open design. The waist belt and hip support structure work together to reduce pressure on the hip region. It improves fit around the waist and distributes force evenly to enhance wearing comfort.

Intended Purpose: Hip orthosis (including abduction orthosis) for external fixation and support of the hip joint.

Indications

- Post-operative use following minimally invasive hip procedures
- Labral repair, with or without gluteus medius repair
- Range of motion (ROM) control following primary or revision hip surgery
- Post-operative use following proximal hamstring repair

Using Instructions

Step 1: Take the left side as an example. Loosen waist belt (#3), thigh straps (#1 & #2), and pad (#5). Thread the left waist belt through the plastic bracket slot (Figures 1a–3a).

Step 2: Connect left and right waist belts through the bracket slot. Adjust to fit waist size. Secure Velcro fastening (Figure 4a).

Step 3: Place pad (#5) in the center. Attach it to Velcro on both sides. Trim back plate (#4) if it interferes (Figure 5a).

Step 4: Wear the brace. Place the frame on the outside of the leg. Align hinge center with the hip joint. Tighten waist belt (Figures 6a–7a).

Step 5: Use pull ring (#6) to adjust waist tightness (Figure 7a).
Step 6: Press button (#7). Adjust connecting rod length. Release to lock (Figures 8a–9a).

Step 7: Tighten thigh straps (#1 and #2). Hold frame during adjustment. Secure all buckles (#10) (Figures 10a–12a).

Step 8: Adjust ROM. Press angle button. Set to required position. Lock buckle (#9) at 90° position. Remove protective pad film after fitting (Figure 10a).

Important Notes

- The connecting rod length is pre-adjusted by the physician to ensure optimal hip stability.
- The pivot center of the angle adjustment mechanism must align with the hip joint center.
- The thigh strap (#1) should be positioned approximately 2 cm below the groin.
- The thigh strap (#2) should be positioned approximately 7 cm above the knee.

Strap Length Adjustment
After determining the required length, wind excess strap material around the strap storage device.

Strap Cutting Adjustment
Determine the required length, open the Velcro at the end of the strap, separate the hook-and-loop fastener, and cut the strap to the desired length. Then reattach the Velcro securely to the strap end.

Abduction / Adduction Adjustment
If required, the clinician may adjust the abduction/adduction angle using the supplied wrench tool.

The patient should follow the physician's instructions regarding hinge angle and ROM settings until medical clearance is given for removal.

Maintenance

Hand wash with mild detergent
Water temperature below 30°C
Air dry only
Do not use dryer or heat
Do not use bleach or softener
Rinse after sweat, sand, or salt exposure

Warning

If you experience pain, swelling, skin irritation, or discomfort while using this product, discontinue use immediately and consult your physician.
This product is a support device only. It is not intended to prevent injury recurrence and does not replace proper medical diagnosis or treatment.

Product materials

- Strap: Nylon
- Bracket buckle: Polyethylene
- Velcro: Nylon
- Connecting rod: Stainless steel

Contraindication

Do not use this product if you have:
Allergy to materials
Skin infection or open wounds
Swelling or ulcers
Severe heart, lung, or bone disease without doctor approval

Storage, Transport, and Shelf Life
1. Shelf life: 5 years
2. Storage conditions: Store in a dry, clean, and well-ventilated environment prior to use.
3. Transport conditions: The product must be protected from dust, moisture, and contamination during storage and transportation

Disposal

Dispose according to local regulations.

Liability

Contact the manufacturer in case of quality problems with the product. Serious accidents that occur during the use of the product need to be reported to the relevant authorities of the EU Member State where the user is located and notified to the manufacturer. The criteria for serious accidents are set out in Article 2, paragraph 65 of the EU Regulation (2017/745).

Warranty

Simple Medical (Shenzhen) Limited warrants to the user who originally purchases this product that it is free from defects in material and workmanship. The sole obligation of Simple Medical (Shenzhen) Limited in the event of breach of warranty shall be to repair or replace the defective product or part(s).

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Simple Medical (Shenzhen) Limited zaručuje původnímu kupujícímu, že tento výrobek je bez vad materiálu a zpracování.